

PALL Life Sciences

Integration of Ion Exchange Chromatography in Downstream Processing of Proteins

SCI Training Course, Cambridge, September 2012

Sylvio Bengio, Ph.D.
Pall Life Sciences

1

PALL Life Sciences

Disclaimer

This presentation is the Confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer.

© 2012 Pall Corporation.

2

PALL Life Sciences

Pall Chromatography Toolbox

Lab Scale Methods & Screening	Process Development	Scale-Up & Production
<ul style="list-style-type: none"> AcroPrep™ and AcroWell™ filter plates. SELDI Services for sorbent screening 1 mL prepacked PRC columns Prepacked AcroSep™ columns Mustang® XT Acrodisc or XTS 	<ul style="list-style-type: none"> Chromatography sorbents Mustang XT5 capsules Mustang XT Acrodisc 5 mL prepacked PRC columns Empty laboratory LRC glass 	<ul style="list-style-type: none"> Chromatography sorbents Mustang XT140 and XT5000 capsules Resolute® columns / Packing support Packing stations PKP and PK systems

3


Defining Ion Exchange Chromatography

▶ Ion exchange chromatography (I-än eks-chänj krö-mä-tä-græ-fë):


“ A separation technique utilizing a stationary phase that contains either acidic groups for exchanging cations or basic groups for exchanging anions present in the mobile phase under defined conditions of pH and ionic strength ”

4

Proteins have different isoelectric points




Isoelectric point (pI)
 Calculated using pK values of amino acids as described in:
 Bjellqvist, et al, *The focusing positions of polypeptides in immobilized pH gradients can be predicted from their amino acid sequences.*
 Electrophoresis 1993, 14, 1023-1031




5


Proteins have different hydrophobicities



GRAVY (Grand Average of Hydrophobicity) score
 Calculated as the sum of hydrophobicity values of all amino acids, divided by the number of residues in the sequence.



6



Ion-Exchange Chromatography

- Routinely used in process-scale protein purification
- Large number of products on the market
- Anion and Cation exchangers
- Applications : Mabs, polyclonal IgG, Recombinant proteins, plasma derivatives, Vaccines ...

8

Main Ligands Used for Ion Exchange


- Anion exchange
 - ♦ Strong **R-N⁺(CH₃)₃**
Quaternary ammonium (Q)
 - ♦ Weak **R-CH₂CH₂NH⁺(CH₂CH₃)₂**
Diethylaminoethyl (DEAE)
- Cation exchange
 - ♦ Strong **R-CH₂SO₃⁻**
Sulfonic acid (S)
 - ♦ Weak **R-CH₂COO⁻**
Carboxymethyl (CM or C)

9


PALL Life Sciences

Formats: Beads or Membrane adsorbers

Sorbents



Membranes



10

PALL Life Sciences


« Salt Tolerant » Ion Exchange Ligands

- Vast majority of proteins acidic. Will be negatively charged around neutral pH or higher, anion exchange (Q, Deae) is more commonly used than cation (S, CM) for capture steps.
- Novel sorbents and membranes allow protein capture at moderate to high conductivities (e.g. >10 mS/cm) with high binding capacities.
 - ◆ For capture : direct load of feedstock
 - ◆ For polishing (flowthrough mode)
- Usually based on primary amine ligands


11

PALL Life Sciences


Typical Steps in Ion Exchange




Equilibration




Load Feedstock




Wash Column



Re-equilibration




Clean in Place





Elution

Confidential

12







PALL Life Sciences

Bind/Elute and Flowthrough Modes

```

graph TD
    A[Target & Contaminants] --> B[Bind and Elute  
Target binds  
Contaminants pass]
    A --> C[Flowthrough Mode  
Target passes  
Contaminants bind]
    
```

16

PALL Life Sciences

Bind/Elute or Flowthrough mode ?

- Relative concentration of target
- Relative pl of target vs main impurity
- Feedstock volume
- Subsequent downstream process
- Target stability to process conditions
- Production facilities and equipment

17

PALL Life Sciences

Selectivity and Resolution

Selectivity (α) : ability of the matrix to separate peaks, i.e. the distance between two peaks. Selectivity can vary from one sorbent to another for the same chemistry

Expressed as : $\alpha = \frac{V_2 - V_0}{V_1 - V_0}$ with V_0 : void volume of the column
 V_1 : elution volume of peak 1
 V_2 : elution volume of peak 2

Resolution : distance between peak maxima compared with the average base width of the 2 peaks. The further apart and the sharper the better the resolution.

Expressed as : $Rs = \frac{2(V_2 - V_1)}{W_1 + W_2}$ with W_1 : peak 1 width at half peak height
 W_2 : peak 2 width at half peak height

18

PALL Life Sciences

Chemical Composition of Sorbents

Polysaccharides Sepharose®, Capto™, Ultrogel®
Sephadex®
HyperCel™, Cellufine™

Synthetic polymers Acrylic
Methacrylate
Polystyrene-DVB

Mineral oxides Silica
Ceramic
Zirconium oxide
Hydroxyapatite

Composite

Trisacryl®, UNOsphere™
Toyopearl® (GigaCap), Fractogel® EMD, Macro-Prep®
Mono®, Source™, Poros®

Spherodex®, Spherosil®
Ceramic HyperD®
HyperZ®, ZirChrom®
CHT™, HA Ultrogel®

PALL Life Sciences

Sorbent Chemistry has an impact on Selectivity

Anion exchangers:

- Q rigid cellulose sorbent
- Rigid Q agarose sorbent
- Q polymeric sorbent

Cation exchangers:

- S rigid cellulose sorbent
- Rigid S agarose sorbent
- S polymeric sorbent

Protein mix eluted by Linear salt gradient

PALL Life Sciences

Capture Step

Media

Media Filtration

Bioreactor

Cell Separation & Clarification

UF/DF

Chromatography Intermediate

Membrane Chrom: DNA Polishing

Final Filtration

Virus Filtration

Final Filtration

Chromatography Capture

21

PALL Life Sciences

Key parameters for Capture

- Selectivity for the target protein
- Dynamic Binding Capacity for the target
- Volume and conductivity of feedstock
 - ◆ Most ion exchangers loose capacity if crude (non-diluted) feedstock is loaded: need for additional unit operations: UF/DF, dilution
- Final column volume and equipment issues

22

PALL Life Sciences

Enhancing Target Protein Capture by IEX


- Check selectivity for your target
- Check « real-life » Capacity for your target
 - ◆ Vendor's specs are only given for guidance, « real-life » capacities tend to be lower !
- Optimize loading, wash/elution pH ...
- Check elution column volumes
- Optimize residence time (start with 4-5 min, decrease to 2 min. if positive)

23

PALL Life Sciences

Prepacked Columns are convenient tools to check both Selectivity and Capacity parameters

- Packing performances guaranteed by the supplier
- Easy to connect to AKTA™ or other system
- Columns can be connected in series to mimic a bed height typical of a pilot-scale column
- Use typical flow rate of 1 mL/min (300 cm/hr)




24

Elution Conditions

- Elution can be achieved as usual by salt gradient.
- Elution may also be done by decreasing the pH. In this case, it allows to recover the product at lower conductivity allowing a direct load to an orthogonal ion exchanger or simplifying formulation

Ion Exchange for Capture : purification of HSA from crude, undiluted plasma



- Direct capture of HSA from undiluted feedstock: Purity > 99%, Yield ~90%
- Elution can be prompted by lowering pH only, no need to add NaCl, this allows a direct orthogonal load on S HyperCel cation exchange sorbent.

High Throughput Sorbent Screening

- HTS methods are well suited to quickly select the best sorbent for a specific step
- Saves time, sample and manpower
- Needed at early stages of process development
- Use 96 well plates combined with Design of Experiments (DoE) software tools



FIRST INTERNATIONAL CONFERENCE
BRNO/CZECH REPUBLIC, OCTOBER 6-7, 2010

Second international conference: Avignon, June 4-7, 2012

High Throughput Sorbent Screening Platform (Pall)




TECAN Freedom EVO[®] workstation allows flexible and precise fully automated high throughput screening of different chromatography sorbents in AcroPrep™ 96-well filter plates.

Capacity of Ion Exchange sorbents


- Dynamic Binding Capacity (DBC) is a key parameter in process chromatography.
 - ♦ To maximize productivity during capture steps, with feedstream having high expression titers (e.g. Monoclonal antibodies 5-10 g/L)
- DBC varies according to the sorbent type, the linear flow rate or residence time
- pH and ionic strength of the feedstock also impact DBC.
 - ♦ New generation of « salt tolerant » sorbents


Dynamic Binding Capacity vs. pH and Conductivity (BSA model protein)




High dynamic binding capacity (DBC) over a wide range of pHs and conductivities at short residence time leads to flexibility and better productivity.

Sorbent: HyperCel STAR AX (Pall) Column: 0.5 cm I.D. x 5 cm bed height (-1 mL); Sample: 5 mg/mL BSA in equilibration buffer; Equilibration buffer: 25 mM Tris-HCl, pH 7.0 – 8.5; Conductivity 3 – 20 mS/cm; Residence time: 1 – 4 min (0.25 – 1 mL/min).








PALL Life Sciences

Ion Exchange at Capture Step: Conclusions

- Check selectivity for the target first
- Choose sorbent that maximizes **Dynamic Binding Capacity** with minimal feed pretreatment (UF/DF, dilution)
- Consider re-equilibration duration, column volumes needed for elution and cleaning in place regimes
- Consider column final volume at production scale and equipment costs

34



PALL Life Sciences

Key parameters for Intermediate step

- Selectivity for the target protein vs main contaminants
 - ♦ e.g. Host Cell Proteins, aggregates, misfolds....
- Consider **Flowthrough mode** (contaminants bind)
- High resolution steps may help in some cases
 - ♦ Smaller beads may also increase process costs

36

PALL Life Sciences

Ion Exchange in Flowthrough mode : Removal of contaminants in rec. Protein (CHO)

Sorbent	DBC for contaminant proteins (mg/mL)
HyperCel STAR AX	~8.5
Q Rigid agarose	~4.5

15 mS/cm

- Target in Flow through, contaminants bind with higher capacity vs conventional sorbent
 - Capacity at 15 mS/cm for contaminants >2-fold higher than Q rigid agarose
- Other data (*not shown*) also demonstrates:
 - Excellent target recovery (no protein loss)
 - Good removal of protein contaminants and pigments
 - Efficient regeneration with 0.1 M HCl followed by 1 M NaOH

CHO Cell Culture Supernatant Courtesy of Cytheris, France

37

PALL Life Sciences

Key parameters for Final Polishing Step


- High Resolution sometimes needed
 - Smaller beads may also increase process costs
- Consider **Membrane adsorbers** instead of packed columns
 - DNA and HCP removal
 - Disposale chromatography


38


PALL Life Sciences

Resolution in IEX is influenced by

- Buffer pH/Protein pl
- Particle Size**
- Load
- Gradient Slope
- Flow Rate







PALL Life Sciences

Membrane versus Resin Chromatography

Packed 40-50 µm Beads
300-1000 Å pores

Mustang™ Membrane
6000 Å pores

Detail
25x 25x
10000 Å = 1 µm

Active surface (green) on polymer matrix (white)
Convective channels Diffusion channels

43

PALL Life Sciences

Benefits of Membrane Adsorbers vs Column Chromatography

- Disposable, plug and play, no time consuming equipment prep
- Less surface area but all of it is accessible to the product
- Orders of magnitude faster flowrates
- Use of a membrane adsorber will allow the polishing step to be accomplished using a device of much reduced volume.
- Ideally suited for capture of large proteins, or polishing

44

PALL Life Sciences

Ion Exchange Mustang® Membranes

- Polyethersulphone (PES) membrane substrate
- Cross-linked polymers containing functional groups – Q, S, E
- Nodular surface - high surface area
- High voids - 0.8 µm pores
- Multiple layers (16) - controlled bed depth


Mustang Q XT Capsules

Mustang Q membrane

45

PALL Life Sciences


Mustang® Membrane Chromatography



Single-Use Mustang Capsules : Q, S

XT Acrodisc®	CLM05	CL3	NP6	NP7	NP8
0.86 mL	10 mL	60 mL	260 mL	520 mL	780 mL

(standard) Acrodisc® Q & S : 0.18 mL, not for scale-up
Mustang E has different volumes : 0.12mL; 10; 40; 160; 320; 480 mL



Re-usable Mustang XT for Q chemistry

XT Acrodisc®	XT5	XT140	XT5000
0.86 mL	5 mL	140 mL	5000 mL

46

ProBioGen
Supporting Biopharmaceutical Growth

Implementation of membrane chromatography as polishing step in IgG1 production

Stefan Franke, Stefan Hartmann, Martin Suhr (ProBioGen AG)

Case study ProBioGen

ProBioGen
Supporting Biopharmaceutical Growth

Introduction

- The presentation focuses on the downstream purification of a monoclonal IgG1 produced in CHO cell culture. Aim of the process development was to establish a GMP-compliant manufacturing process at 250 L scale culture for supply of material for toxicological studies and later-on for phase I clinical trials.

45

Outline ProBioGen
Supporting Biopharmaceutical Growth

- I. Evaluation of mixed-mode chromatography as a polishing step for the removal of MAb aggregates.
- II. Implementation of membrane chromatography as a polishing step for host cell DNA removal: from development to production.
- III. Removal of HCP during the membrane chromatography step.
- IV. Pilot-scale study of viral clearance.

49

Human IgG1 from CHO culture - DSP Process ProBioGen
Supporting Biopharmaceutical Growth

- Three basic chromatography steps.
- Two virus depletion steps: low pH treatment and nanofiltration.
- Polishing (DNA and HCP removal) - Two techniques were investigated:
 - MEP HyperCel™ mixed-mode chromatography for aggregate removal.
 - Anion exchange Mustang® Q chromatography membrane adsorber (AIEX) used in negative-mode

```


graph TD
    A[Cell free supernatant] --> B[Capture Protein A]
    B --> C[Intermediate CIEX]
    C --> D[Polishing MEP AIEX]
    D --> E[Pure Antibody Batch]
  
```

50

Material and Methods ProBioGen
Supporting Biopharmaceutical Growth

Material:
IgG1 was expressed in serum-free CHO cell culture (fed-batch).

Methods:
All chromatography steps were performed using ÄKTA™ systems (GE HealthCare). MEP HyperCel™ sorbent and Mustang® Q membrane adsorbers were supplied by Pall.




Analytics:

- IgG1 purity: Size Exclusion Chromatography
- Residual DNA: Picogreen Assay (Invitrogen) or quantitative PCR.
- Residual HCP: generic anti CHO-HCP-ELISA (Cygnus Technologies).
- Virus spiking study: external service provider.

51

IgG1 recovery after Mustang Q Polishing



Unit: Scale dependent
 Buffer A: 20 mM Na-Phosphate pH 7.0
 Buffer B (optional for eluting): 20 mM Na-Phosphate pH 7.0, 1 M NaCl
 Conditioning of load: 1:2 dilution of the CIEX eluate with Buffer A


The 10 ml Mustang Q capsule was selected for production at the 50 L scale, and the 60 ml capsule for 250 L.

Run	Process Volume (mL)	Mustang Q unit	Conductivity	IgG1 concentration (UV280nm)	Protein Load (g product/ml bed volume)	Mab recovery
Consistency run 02	1010	10 ml. Capsule	11.3 mS/cm	1.7	n/a	98.5 %
250 L Run	37600	60 ml. Capsule	12.9 mS/cm	1.7	1.07	96.1 %
50 L Run	5609	10 ml. Capsule	11.8 mS/cm	2.4	1.34	97.9 %
250 L GMP run	40150	60 ml. Capsule	12.5 mS/cm	1.7	1.14	104.2 %

Data showed high IgG1 recovery (> 96 %) independent of the scale applied.

55


Process conditions of the membrane adsorber for production



	Tox 1	Tox 2	GMP
Scale	250 L	50 L	250 L
Process volume	225 L	44 L	238 L
Process conditions	Mustang Q capsule: MV: 60 mL Flow rate: 36 L/h (10 MV/min)*	Mustang Q capsule: MV: 10 mL Flow rate: 30 ml/min** (3 MV/min)*	Mustang Q capsule: MV: 60 mL Flow rate: 36 L/h (10 MV/min)*


*MV: membrane volume according supplier specification
 **: due to technical reasons as pressure maximum was reached

Mustang Q 60 mL Capsule



56

Impurities: Removal of residual DNA



Results of Host Cell DNA Depletion (Picogreen Assay)

Process Step	Development Run 1 (10 L Scale)		Development Run 2 (10 L Scale)	
	DNA ng/ml	Product mg/ml	DNA ng/ml	Product mg/ml
Cell culture supernatant	5951	n/a	7538	n/a
After 1st chromatography (Protein A)	19	6.6	22.2	7.2
After 2nd chromatography step (CIEX)	4.2	3.4	2.2	3.4

After CIEX the DNA content was already in the low ng range. The Picogreen Assay could not be applied for analysis of AIEX samples so samples were analyzed externally using a qualified qPCR method (assay sensitivity < 1pg).

Results of Host Cell DNA data determined by qPCR

Sample/Process step	Tox - 250 L	Tox - 50 L	GMP - 250 L
Final bulk	< 40 pg/ml*	< 17 pg/ml*	< 16 pg/ml*

*: below lower limit of quantification of the assay (LLOQ)

The membrane adsorber step showed effective removal of residual DNA at production scale.

57

Conclusions and Summary

- A membrane adsorber step was successfully implemented as single-use polishing step for removal of residual DNA in a IgG1 purification process at production scale.
- Selection of the membrane adsorber format (60 mL Mustang Q capsule for 250L culture) took into account a safety factor (potential biological variation of the feedstock composition).
- Effective MuLV virus removal using Mustang Q adsorber was shown in a scale-down model while MVM virus removal needed further optimization (load conductivity).
- Overall virus clearance capacity was sufficient for manufacturing of a safe drug substance for clinical trials.

58

PALL Life Sciences

Thank You !
sylvio_bengio@europe.pall.com



59
