


## Ion Exchange and Mixed Mode Chromatography in Process Separations

SCI 2012 Technical Training  
Cambridge, September 2012

Dr. Sylvio Bengio



---

---

---

---

---

---

---

---

### Disclaimer

This presentation is the copyright work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer.

© 2012 Pall Corporation.



2

---

---

---

---

---

---

---

---

### Outline

- Mixed Mode Chromatography : principles
- Purification of Monoclonal Antibody using a Two-Step Scheme without Protein-A
- Orthogonal separations: mixed mode and ion exchange chromatography, two case studies



3

---

---

---

---

---

---

---

---

## Mixed-Mode Chromatography

- Exploits multiple, distinct protein-ligand interactions to adsorb target proteins or impurities.
- Offers new solutions where traditional chromatographic methods are not effective.
  - Where feedstream conductivity is too high for efficient capture on traditional ion exchange resins
  - Alternative/complement to conventional Hydrophobic Interaction (HIC) or hydroxyapatite
  - Separations where affinity ligands are too expensive

---

---

---

---


---

---

---

---

## Mikhail Semenovich Tswett (1872-1919) the « inventor » of chromatography



" Like light rays in the spectrum, the different components of a pigment mixture, obeying a law, are resolved on the calcium carbonate column and then can be qualitatively and quantitatively determined. I call such a preparation a chromatogram, and the corresponding method the chromatographic method."

---

---

---

---

---

---

---

---

## Mixed-Mode Chromatography

### Mixed Mode Chromatography A Multimodal Separation Technique for Biopharmaceutical Purification

Wednesday 27 June 2007  
SCI, London, UK



---

---

---

---

---

---

---

---

### Commercial Mixed Mode or Multi-Mode Ligands

- Hydroxyapatite (calcium phosphate)
- Trichlorotriazine dyes
- Multi-modal cation exchange (N-Benzyl-homocysteine)
- Multi-modal anion exchange (N-Benzyl-N methyl ethanolamine)
- 4-MEP (4-Mercapto Ethyl Pyridine)
- HEA (Hexylamine)
- PPA (Phenyl Propylamine)

7

---

---

---

---

---

---

---

---

### Key Sorbent design Parameters

**Ligand density:** Cooperativity between ligands for protein adsorption in physiological conditions dictates high densities.

**Hydrophobicity:** Since this parameter mediates the adsorption, the nature of the linker or of the activating agent is critical.

**pK of the ligand:** This parameter influences the pH of elution of the antibody. To have mild desorption conditions this value has to be between 4 and 9.

8

---

---

---

---

---

---

---

---

### MEP, HEA, PPA HyperCel™

- Binding:** dominant mode is hydrophobic interaction.
  - Binding typically achieved without addition of binding-promoting salt, or at significantly lower salt concentration than for conventional HIC.
- Desorption:** driven by electrostatic charge repulsion.
  - Target pool typically recovered in dilute or relatively low conductivity buffer.

MEP

\*OCCNCCSCC[n+]1ccccc1

HEA

\*CCN

PPA

\*CCCNc1ccccc1

9

---

---

---

---

---

---

---

---

**PALL Life Sciences**

### Interaction of 4-MEP Ligand with Antibody or Protein

**Adsorption at near-neutral pH, physiological conditions**

Hydrophobic Interaction

$pK_a = 4.8$

or = non-antibody = protein

**Desorption at pH 4.0 - 5.8**

Electrostatic repulsion

pH	% in (+) form
4.8	50%
5.8	10%

10

---

---

---

---

---

---

---

---

**PALL Life Sciences**

### Proteins have different charges and hydrophobicities

**Isoelectric point**

**GRAVY score**

11

---

---

---

---

---

---

---

---

**PALL Life Sciences**

### Optimization of Elution-pH on Mixed-Mode Sorbents

- pH 5.5
- pH 5.2
- pH 4.9
- pH 4.6
- pH 4.3
- pH 4.0
- pH 3.0

**Elution of relatively basic proteins.**

**Elution of relatively hydrophilic proteins.**

**Elution of relatively acidic proteins.**

**Elution of relatively hydrophobic proteins.**

**Elution of basic / hydrophilic impurities**

**Elution of target antibody or other target protein**

**Elution of acidic / hydrophobic impurities**

12

---

---

---


---

---

---

---

---




---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---

### Mixed-Mode Chromatography for Post-Capture Purification of Monoclonal Antibodies

- Capture on Protein A sorbent as a standard method ; Cation exchange in some cases
- Post-capture typically done by IEX or HIC
- Mixed-mode sorbents of growing interest because

  - Orthogonal to other methods (cation exchange, hydroxyapatite)
  - Optimization allows selective desorption of MAb from contaminants (HCPs, aggregates).
  - Contributes to decrease purification costs.

**PALL** Life Sciences

---

---

---

---

---

---

---

---

## Scale-Up of Mixed Mode Chromatography



Scale-up on MEP HyperCel sorbent, presented by Dr. Richard Francis, SCI Mixed Mode Conference, London, June 2007


 © Pall Corporation 2012

---

---

---

---

---


---

---

---

## Industry Case Studies

- 1. Two-Column Processes to Purify Monoclonal Antibodies Without use of Protein-A Sorbent
- 2. Combination of mixed mode chromatography and membrane adsorbers


 © Pall Corporation 2012

---

---

---

---

---

---

---

---

## Two-Column Processes to Purify Monoclonal Antibodies Without use of Protein-A Sorbent

UFDF

Mixed Mode on MEP HyperCel

Viral Inactivation (Low-pH)

Cation Exchange Chromatography

Viral Filtration

• Prep for CIEC  
 • UFDF in both schemes allows comparison of MM & CIEC

UFDF


Cation Exchange Chromatography

Viral Inactivation (Low-pH)

Mixed Mode on MEP HyperCel

Viral Filtration

A Two-Column Process to Purify Antibodies Without Protein-A, G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Figure-3.


 © Pall Corporation

---

---

---

---

---

---

---

---

### Purification of Monoclonal Antibody using a Two-Step Scheme without Protein-A

Mode	In Eluate or Load	HuMAb-1	HuMAb-2
Mixed Mode on MEP HyperCel	HCP in load (ng/mg)	323,279.	157,069.
	HCP in eluate (ng/mg)	3,470.	1,810.
	DNA in load (pg/mg)	138,091.	783.
	DNA in eluate (pg/mg)	4.12	0.19
	Purity (% monomer in eluate)	99.91	98.82
	Step Recovery (%)	83.	84.
Cation Exchange	HCP in eluate (ng/mg)	68.	24.
	DNA in eluate (pg/mg)	0.4	0.2
	Purity (% monomer in eluate)	100.	99.82
	Step Recovery (%)	85.	84.
Overall Recovery (%)		71.	71.

"A Two-Column Process to Purify Antibodies Without Protein-A", G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-1.

© Pall Corporation

PALL Life Sciences

19

### Purification of Monoclonal Antibody using a Two-Step Scheme without Protein-A

Mode	In Eluate or Load	HuMAb-1	HuMAb-2
Mixed Mode on MEP HyperCel	HCP in load (ng/mg)	323,279.	157,069.
	HCP in eluate (ng/mg)	3,470.	1,810.
	DNA in load (pg/mg)	138,091.	783.
	DNA in eluate (pg/mg)	4.12	0.19
	Purity (% monomer in eluate)	99.91	98.82
	Step Recovery (%)	83.	84.
Cation Exchange	HCP in eluate (ng/mg)	68.	24.
	DNA in eluate (pg/mg)	0.4	0.2
	Purity (% monomer in eluate)	100.	99.82
	Step Recovery (%)	85.	84.
Overall Recovery (%)		71.	71.

"A Two-Column Process to Purify Antibodies Without Protein-A", G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-1.

© Pall Corporation

PALL Life Sciences

20

### Purification of Monoclonal Antibody using a Two-Step Scheme without Protein-A

Mode	In Eluate or Load	HuMAb-1	HuMAb-2
Mixed Mode on MEP HyperCel	HCP in load (ng/mg)	323,279.	157,069.
	HCP in eluate (ng/mg)	3,470.	1,810.
	DNA in load (pg/mg)	138,091.	783.
	DNA in eluate (pg/mg)	4.12	0.19
	Purity (% monomer in eluate)	99.91	98.82
	Step Recovery (%)	83.	84.
Cation Exchange	HCP in eluate (ng/mg)	68.	24.
	DNA in eluate (pg/mg)	0.4	0.2
	Purity (% monomer in eluate)	100.	99.82
	Step Recovery (%)	85.	84.
Overall Recovery (%)		71.	71.

"A Two-Column Process to Purify Antibodies Without Protein-A", G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-1.

© Pall Corporation

PALL Life Sciences

21

"A Two-Column Process to Purify Antibodies Without Protein-A", G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-2.

Mode	In Eluate or Load	HuMab-6	HuMab-7
Cation Exchange	Column Volume (L)	5.3	5.0
	HCP in load (ng/mg)	54,502.	251,467.
	HCP in eluate (ng/mg)	1,047.	379.
	DNA in load (pg/mg)	nd	1,140,000.
	DNA in eluate (pg/mg)	nd	20
	Purity (% monomer in eluate)	98.47	99.87
	Step Recovery (%)	84.	88.
Mixed Mode on MEP HyperCel	Column Volume (L)	10.6	4.7
	HCP in load (ng/mg)	1047.	366.
	HCP in eluate (ng/mg)	39.	11.
	DNA in load (pg/mg)	nd	20.
	DNA in eluate (pg/mg)	< 6	< 1
	Purity (% monomer in eluate)	98.85	99.94
	Step Recovery (%)	95.	92.
Overall Recovery (%)		80.	81.

© Bell Corporation



"A Two-Column Process to Purify Antibodies Without Protein-A", G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-2

Mode	In Eluate or Load	HuMab-6	HuMab-7
Cation Exchange	Column Volume (L)	5.3	5.0
	HCP in load (ng/mg)	54,502.	251,467.
	HCP in eluate (ng/mg)	1,047.	379.
	DNA in load (pg/mg)	nd	1,140,000.
	DNA in eluate (pg/mg)	nd	20
	Purity (% monomer in eluate)	98.47	99.87
	Step Recovery (%)	84.	88.
Mixed Mode on MEP HyperCel	Column Volume (L)	10.6	4.7
	HCP in load (ng/mg)	1047.	366.
	HCP in eluate (ng/mg)	39.	11.
	DNA in load (pg/mg)	nd	20
	DNA in eluate (pg/mg)	< 6	< 1.
	Purity (% monomer in eluate)	98.85	99.94
	Step Recovery (%)	95.	92.
Overall Recovery (%)		80.	81.

© Pall Corporation



**"A Two-Column Process to Purify Antibodies Without Protein-A", G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-2**

Mode	In Eluate or Load	HuMab-6	HuMab-7
Cation Exchange	Column Volume (L)	5.3	5.0
	HCP in load (ng/mg)	54,502.	251,467.
	HCP in eluate (ng/mg)	1,047.	379.
	DNA in load (pg/mg)	nd	1,440,000.
	DNA in eluate (pg/mg)	nd	20
	Purity (% monomer in eluate)	98.47	99.87
	Step Recovery (%)	84.	88.
Mixed Mode on MEP HyperCel	Column Volume (L)	10.6	4.7
	HCP in load (ng/mg)	1047.	366.
	HCP in eluate (ng/mg)	39.	11.
	DNA in load (pg/mg)	nd	20.
	DNA in eluate (pg/mg)	< 6	< 1
	Purity (% monomer in eluate)	98.85	99.94
	Step Recovery (%)	95.	92.
Overall Recovery (%)		80.	81.

© Pall Corporation





"A Two-Column Process to Purify Antibodies Without Protein-A", G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-2.

© Bell Corporation

 Life Sciences

	HuMAB-6 Log <sub>10</sub> PFU Reduction	HuMAB-7 Log <sub>10</sub> PFU Reduction
Viral Challenge	<ul style="list-style-type: none"> <li>• 2000 mg dose</li> <li>• 9.87 logs infectious &amp; non-infectious viral particles</li> </ul>	<ul style="list-style-type: none"> <li>• 1000 mg dose</li> <li>• 9.33 logs infectious &amp; non-infectious viral particles</li> </ul>
Cation Exchange Chromatography on Fractogel SE Hicap	2.23	2.95
Low pH Treatment	5.04	>5.69
Mixed Mode Chromatography on MEP HyperCel	5.15	>4.59
Viral Filtration	>4.92	>5.95
<b>TOTAL REDUCTION</b>	<b>&gt;17.34</b>	<b>&gt;19.18</b>

A Two-Column Process to Purify Antibodies Without Protein-A, G.M. Ferreira, J. Dembecki, A. Patel, A. Arunakumari, BioPharm International, May 1, 2007. See Table-3.

**PALL** Life Sciences

## PX'Therapeutics Presentation

### Case Study:

### Optimisation of mixed-mode intermediate purification on *E. coli* protein and subsequent membrane polishing step

N. Maguet

Presented at SBCN 2010, october, 19<sup>th</sup>  
Ecole Centrale de Lyon, Ecully, France

**Featuring MEP HyperCel™ and Mustang® Q**

Courtesy of Nicolas Maguet, PX'Therapeutics  
France



Disclaimer. This presentation is the Confidential work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer. © 2011 Pall Corporation



---

---

---

---

---

---

---

---

---

---

---

---

---

---


---

---

---

---

[sylvio\\_bengio@europe.pall.com](mailto:sylvio_bengio@europe.pall.com)



---

---

---

---

---

---