



Bloom Growth in Chocolate: The Application of Profilometry to Study Transformations and Aging

Paul Smith

Cargill R&D Centre Europe, Vilvoorde, Belgium

Dérick Rousseau

Ryerson University, Toronto, Canada

Application of Profilometry

- **Chocolate Bloom**
- **Other techniques**
 - Atomic Force Microscopy
 - Environmental Scanning Electron Microscopy
- **Profilometry**
- **Aging of Chocolate**
- **Conclusions**

Chocolate Bloom



- **Loss of gloss and graying of chocolate surface**
- **Poor appearance – Any reason**
- **2 main causes**
 - Sugar bloom
 - Sugar recrystallisation, humidity problems
 - Fat bloom
 - Cocoa butter recrystallisation
 - Transformation from desired form V to stable form VI
 - Poor tempering
 - Poor storage
 - Filling fats

Chocolate Bloom Studies


- **Difficult to study and Monitor Development**
- **Surface behaviour**
- **Nature of chocolate leads to difficulties for microscopy**
- **Slow time development**
- **Difficult to observe changes**

- **Development of Different techniques**

Principle of AFM




Fresh Dark Chocolate




06171642.002

Fresh Dark Chocolate (3D view)

Digital Instruments NanoScope
Scan size 20.00 μm
Scan rate 0.7825 Hz
Number of samples 512
Image Data Height
Data scale 800.0 nm



□ view angle
☀ light angle




X 5.000 $\mu\text{m}/\text{div}$
Z 800.000 nm/div


06171642.002

Dark Chocolate after 2 months

Digital Instruments NanoScope
Scan size 20.00 μm
Scan rate 1.001 Hz
Number of samples 512
Image Data Height
Data scale 3.000 μm



Dark Chocolate after 5 Months (Bloom Crystal)




Atomic Force Microscopy Conclusions


- **Resolution of very small structures**
- **Physical information about surface**
- **No special treatment**

- **Contact technique**
- **Destructive**
- **Time consuming**

SEM v ESEM




SEM




ESEM

Dark Chocolate Surface




Cross Section of Chocolate

Bubbles



Bloomed Chocolate



Environmental Scanning Electron Microscopy

Conclusions

- **Very good resolution possible**
- **Surfaces do not need to be as well defined as for AFM**
- **Cross sections can be measured**
- **Whole samples can be studied and repeat measurement is possible.**

- **Samples need to be placed in an environmental chamber**
- **Risk of damage by electrons**

Profilometer

Optical
Microscope



Computer


Profilometer: Principles

- **Interference patterns are measured at different heights through the range of the surface.**
- **From these the computer calculates a 3D surface map.**
- **Stitching of images can give views of very large surface areas.**


Melting of Chocolate

- Heat at 1°C per minute from 20°C
- Live measurement


Milk chocolate 0 minutes 20°C




Milk chocolate 5 minutes 25°C




Milk chocolate 7minutes 27°C




Milk chocolate 10 minutes 30°C



Milk Chocolate 13 minutes 33°C




Milk Chocolate 15 minutes 35°C




Chocolate Storage

- **Stored over 2 weeks**
- **Isothermal at 20°C**
- **Temperature cycling 6 hours at 30°C & 6 hours at 20°C**


Milk chocolate 0 days




Milk Chocolate 0 days




Isothermal Chocolate 13 days




Isothermal Chocolate 13 days



Cycled Chocolate 13 days



Cycled Chocolate 13 days



Roughness Changes of the Samples

Sample	Rq	Z range (nm)
Day 0	119.80	3.43
Day 13 Isothermal	152.62	3.83
Day 13 Cycled	482.53	5.41

Profilometry Conclusions

- **Good surface representation is possible**
- **No strange treatments are needed**
- **The same spot can be studied time after time**
- **Quick and easy technique**

- **High quality surfaces are needed for good resolution**

Conclusions

- **All 3 techniques are useful to understand the structure of chocolate and the nature of chocolate bloom**
- **In particular profilometry is useful for studying surface changes and monitoring transformations in real time**

Acknowledgements

- **Annika Dahlman & YKI Institute for Surface Chemistry, Stockholm, Sweden.**